

# झापा गाउँपालिका गाउँ कार्यपालिकाको कार्यालय



प.सं.:२०६१/६२ च.नं: 2629 टाघनडुब्बा झाण, कोशी प्रदेश, नेपाल Email: jhapagaunpalika@gmali.com | website: jhapamun.gov.np

> मितिः२०८९/१२/१० ने.सं-११४५, चिल्लागा, ९ आइतबार

श्री इजाजत प्राप्त फर्महरू सबै

विषयः दररेट पेश गर्ने सम्बन्धमा ।

प्रस्तुत विषयमा अपाङ्गता भएका व्यक्तिका लागि आवश्यक संलग्न स्पेशिफिकेशन बमोजिमका सहयोगी सामग्री के कित दरमा उपलब्ध गराउन सक्नुहुन्छ? सोको दररेट ७ दिनभित्र पेश गर्नुहुन अनुरोध छ ।

गजेन्द्र थपलिया

प्रमुख प्रशासकीय अधिकृत

गजेन्द्र थयलिया प्रमुख प्रशासकीय अधिकृत

मितिः २०८१/१२/०५




ASSISTIVE DEVICE	SPECIFICATION	Quantity	Rate	REMARKS
ASSISTIVE DEVICE	Amplification frequency response and maximum output (OSPL) that meets requirements of standard prescription formul (e.g. Desired Sensation Level v.5 or National Acoustics Laboratories NL2)	ae		Audio service Germany
Semi-Digital and Digital hearing aids with volume control	Frequency range should be 200–4500 Hz (minimum)  Digital Hearing aids with programs, memoral control etc.  Suitable for any severity of hearing loss from moderate (40 dB) to profound (90 dB) degree			Starkey, American
	At least two preset hearing programmes Preconfigured (stock or disposable) ear inse in variety of sizes or custom made earmould	rts s		
2 Walker	Rigid steel frame Foldable with net height of cm Foam hand grips for extracomfort. Width/ folded depth: 52cm/10cm Height: From 79 to 91.5 Weight bearing capacity: 3.3kg/125kg	a 1		
3 Cane/ L-stick	Size: Large, Medium, Small  Constructed from alumin alloy tubes.  ABS handle with comfor grip.  Max/min adjustable heig 780/742 ± 5 mm  Height-adjustment mechanism: Telescopic ,quic release  Weight bearing capacity 110kg  Size: Large, Medium, Sm	table ht:  k- : all		
4 Tripod stick	Constructed from mild saluminium alloy tubes. Anodized for long lastin broader base Soft gel on handle top Tips securely fitted, non and replaceable or interchang Max/min adjustable heig 984/745 ± 5 mm Height-adjustment mechanism: telescopic, quick release Weight bearing capacity 110kg Size: Large, Medium, Sr	g and -slip, geable ght:		A STATE OF THE STA

Axillary crutch		Constructed from mild steel aluminium alloy tubes. Double molded hand grip with soft cushioning.  Height-adjustment mechanism: double height adjustment (handgrip and shaft)  Tips securely fitted and made from durable rubber  Max/min adjustable height: 984/745 ± 5 mm  Weight bearing capacity: 120kg  Size: Large, Medium, Small		Charles of the control of the contro
6 Elbow crutch	Flipping should be should	Shaft: • height-adjustable • fitted with tip  Height-adjustment mechanism: • quick-release • single height adjustment (shaft only) • double height adjustment (cuff and shaft)	1	
Quadripod		Bariatric range for users over 120 kg  Constructed from steel tubes.  Anodized for long lasting and broader four leg base Soft gel on handle top Tips securely fitted, non-slip, and replaceable or interchangeable Max/min adjustable height: 984/745 ± 5 mm Height-adjustment mechanism: telescopic, quick-release Weight bearing capacity: 110kg Size: Large, Medium, Small	1	
		Frame Material MS Chrome  Seat Width 18.5 Inches  17.5 - 21.5 Inches  Height Yes  Adjustable Model Weight 7 KG		JUST STORM

	Frame Style Foldable User Capacity 100 KG Wheel No  Fixed or removable backrest Seat with aperture Removeable collection bucket (pan)	के जिस गाउँपार राष्ट्रमालकाको क के जुक्का स्ट्राप
Wheel chair with cushion  Semi-active wheelchair	Durable metal body Foldable Flip up arm rest Swing away and detachable foot rest Quick release rear wheels Solid PU Tyer Fixed back Wheelchair seat cushion	200 at 100 at 10
	Cushion design:  overall height approximately 70–120 mm  allows full or deep offloading or immersion (sinking into) and envelopment (conforming to shape) 40–45 mm) or deeper	FREEDOM 5000
	Cushion size range:  • should match wheelchairs provided  • seat width range in 25–50 mm increment options or adjustable  • different length options or adjustable in length	
	Supplied with cover that meets general design requirements; for users with incontinence, a water-resistant cover or lining to protect cushion should be supplied  Typical user	
	<ul> <li>Child or adult with intermediate or advanced posture support needs who requires variable tilt or recline, or larger range of static tilt or recline than available on manual wheelchairs</li> </ul>	
	<ul> <li>For users who self-propel or need to be assisted</li> <li>For use in both indoor and uneven urban, peri-urban and rural outdoor environments</li> <li>Specific characteristics</li> </ul>	
	Wheelchair with folding or rigid frame, three or four wheels, push-handles, rigid seat and backrest, armrests or tray table, and footrests     Range of integrated or included PSDs and posture support cushion	
	Large-range adjustable recline and large-range adjustable and variable (dynamic) tilt      Rear wheels and front castors appropriate for indoor and outdoor mixed terrain	
	For intermediate and advanced service levels only: Frame:  Frame with mechanism(s) to fold or	

 Rigid backrest and seat; removable without tools on folding frame

## Requirements for standard configuration



- Adjustable recline18 with minimum 20 ° with tools; optional variable (adjustable without tools) recline minimum 20 °
- integrated mechanism or adjustment options to prevent tip-over when chair used with tilt or recline
- Footrests with adjustable height; available in at least two adjustment ranges on all sizes of wheelchair to accommodate people with shorter and longer legs; two flip-up, swing-away or removable footrests for wheelchairs with folding
- Mechanism or adjustment options available to maintain backwards stability with maximum combined seat and backrest configuration adjustments

for self-propelled wheelchairs: adjustment options to improve reach to rear

### Rear wheels:

10 Wheelchair with postural support (CP wheelchair)

- diameter, appropriate for size of wheelchair, e.g. 559–660 mm for self- propelled wheelchairs or from 305 mm for assistant-propelled wheelchairs
- quick-release or removable without tools on rigid frames for self-propelled wheelchairs

#### Castors:

- diameter e.g. 76-203 mm
- width e.g. 13-50 mm
- may be larger and wider if only one castor (depending on terrain) Range of size and width options available
- Puncture-proof or pneumatic rear tyres and front castors

#### Back rest

- height adjustable or available in different lengths to optimize trunk support (at least at shoulder height or higher)
- not integrated into back-posts; separate component removable from back-posts if required, or part of integrated, removable seat—backrest unit

े गाउँपारिका विभागितकाकी विभागितकाकी विभागितकाकी विभागितकाकी विभागितकाकी विभागितकाकी विभागितकाकी विभागितकाकी विभागितकाकी विभागितकाकी

1

A STUDIE STEEDE

	optional: rigid backrest independently adjustable from back-posts, with forward and backward and angle adjustments can be mounted at different heights on back posts		के किया जा के किया के किया के किया किया किया किया किया किया किया किया
	<ul> <li>headrest attached to back of backrest and with adjustable forward and backward, sideways, height and angle; must be able to support head behind, in line or in front of back support as required</li> </ul>		केरेशी प्रदेख
	<ul> <li>pelvis strap (positioning belt) with adjustable length</li> <li>adjustable and removable knee separator pad (medial knee support); may be integrated in cushion; standard on child size and optional on</li> </ul>		
	adult size  calf or foot straps with adjustable length  tray table or flip-up or removable armrests with adjustable height		
	outside (lateral) thigh or knee pads (supports) for seat, shoulder harness and chest strap		
	Externally applied customized appliances/ device used to support, assist, or correct the weakened segment of the lower extremity in the subtalar region, ankle, and foot.		
11 Ankle-Foot Orthotics (AFO) / Orthosis for foot drop	Specification  1. Made of durable and elastic polypropylene.  2. Suitable for all types of shoes, can be used	1	
	while wearing shoes  3. Help with dorsal flexion  4. Light weight, strong durability and load bearing  5. Colors and details can be customized		
	Size: Individual measurement should be taken		
	Includes: CTEV shoes, Shoe raise		
	Orthopedic shoes are specially designed footwear intended to provide support, comfort, and relief for individuals with various foot and lower limb conditions or deformities. These shoes are designed to address specific orthopedic issues, promote proper foot alignment, and enhance overall foot health.		
	Orthopedic shoes come in various styles and can accommodate a range of foot conditions, from common problems like bunions and hammertoes to		Sett at
12 Modified Shoes	more serious conditions associated with arthritis or diabetes	1	STORE RELEGIES STREET

	Specification  Materials most include leather, foam rubbers, cellular polymers having long durability and comfortable wearing with ease.  Size: Individual measurement should be taken	े जा
	CP Chair (with adjustable incline &tray) Hardwood & ply board constructed, naturally finished Chair reclines through an arc and can be adjusted to any desired position.  Comes with adjustable foam	
3 CP Chair	Padded headrest which can be formed & reformed as needed.  High back supports & holds the child with wide belt.  Laminated activity Tray size 56cm x 63cm with polished wooden raised edges, holds the items from falling.	
	Adjustable Seat &Foot board are 30cm wide and suits the small children. Mounted on ball bearing castors.      Knee separator  Size: Individual measurement should be taken	
	Weight: 35kg  Material: PVC+ Wood and steel frame+ Soft cushion + Vinyl nylon straps  Suitable age: 5-15 years old children (90cm - 130cm)  Use scenario: Applicable in classrooms, rehabilitation centers, therapy rooms at home	
14 CP standing frame	able Height Adjustable for Customized Support     Provides Proper Standing Positions     Improves Posture, Trunk Control, and Balance     Supports Weight-Bearing Exercises and Lower Limb Strengthening     Durable Construction for Longevity	1
	Size: Individual measurement should be taken  A KAFO is a device that is used to control instabilities in the knee and lower limb by maintaining proper alignment and controlling	THE THE LEGIST SHE

	5. Colors and details can be customized		
	Size: Individual measurement should be taken		ीत माउँपा है। कार्यपालिकाकी टायनहरूवा, क्रार्थ
16 Rollator	Short- or long-term use for adult who is frail or needs support to balance or bear weight through their legs Specific characteristics  1. Pushed in front of user 2. Frame: • folding or rigid • height-adjustable handgrips • four wheels • seats 3.Height-adjustment mechanism: • telescopic • quick-release 4.Handgrips: • moulded grips  Requirements for standard configuration • Two handgrips • Four wheels with diameter 100 mm • Driving and parking brakes • Fold-down seat May include basket Size range for users over 120 kg	1	
17 White Cane	Aluminum alloy Made of high quality 4 section aluminum construction. Canes fold and unfold with ease, good to storage in your bag. Cane shaft is covered with red and white reflector tape for nighttime visibility. It comes with standard wrist strap and reinforced nylon tip Weight of 0.53 lb. Length range: Tall adults: 145 cm ± 5 cm Adults: 120–140 cm ± 5 cm Children: 70 cm ± 5 cm Plastic or rubber handgrip of length 20 cm and diameter 2.5 cm and wrist loop  Upper/ Lower limb		
	prosthesis  Size: Individual measurement should be taken		Walter State of the State of th

8 Prosthesis		1	ति माउँपाटि जार्थमातिकाको उत्तर रोधनङ्ख्या, अर्गार्थ रोधनङ्ख्या, अर्गार्थ
19 Magnifying glass	Convex lens that is used to produce a magnified image of an object	1	-70
Air Mattress  20	Effective alternating pressure system and is an ideal for prevention and bedsores treatment of patients at low risk of pressure ulcers.  This alternating pressure pad is composed of 2.5"-high bubble cells overlay mattresses and easy to install.  It brings to you anti- decubitus air pump and bubble mattress specifically designed for prevention of bedsores.  Simple to operate, portable and lightweight, this mattress useful for pressure ulcer prevention and treatment.  On/Off switch & rotary dial for pressure setting.  Comfortable, Perfect finish, Robust design, Stunning design.  AIR PUMP  Power supply:  AC220V/50H  Z  Air output: 6 - 7 liter/min.  Pressure range: 40-100 mmHg  Cycle time: 5 to 6 min.  Power consumption:  26.5cm (L) x  11cm (W)  x 12.5cm (H)  Fuse: 1A  Weight: 1.5 kg  Noise: 45db	1	

y wool

William Wall Britished